

Characterization of biopolymers by advanced polymer chromatography (APC)

Magdolna R. Mihályi
Research Centre for Natural Sciences

Project meeting
„Joint chemical laboratory for the service of bioeconomy in the Slovak-Hungarian border region”
Interreg, SKHU/1902/4.1/001/Bioeconomy

Research Centre for Natural Sciences,
Magyar tudósok körútja 2, Budapest H-1117, Hungary

22 September, 2021



Building Partnership



SEC/GPC/APC

SEC: Size Exclusion Chromatography (1959, Porath and Flodin)

GPC: Gel-Permeation Chromatography (1974, Down Chemical. Co.)

stationary phase: synthetic polymer, e.g. PS

APC: Advanced Polymer Chromatography (2004, Waters Co., UPLC)

stationary phase: rigid, 2.5 μm -size modified silica particles with pore size of 45Å - 900Å.

SEC/GPC/APC

- molar mass averages,
- molar mass distribution of synthetic and biopolymers

number average molecular weight

weight average molecular weight

polydispersity index


$$M_n = \frac{\sum N_i M_i}{\sum N_i},$$

$$M_w = \frac{\sum N_i M_i^2}{\sum N_i M_i}$$

$$D = \frac{M_w}{M_n}$$


SEC principle

- Polymers are separated by hydrodynamic volume
- Big One Comes Out First (BOCOF) followed by the smaller molecules



APC columns for aqueous and organic polymer separation

Ethylene Bridged Hybrid (BEH) technology, Waters



- strong and rigid particles
- particle size: 1.7 and 2.5 μm
- resist shrinking, swelling
- easy solvent switching
- high reproducibility

APC columns

10 small columns, diameter: 4.6 mm; length: 150 mm

	Solvent	Temp. limit (°C)	pH	Pore size (Å)	Particle size (µm)	Linear range (g/mole)
ACQUITY APC XT 45	organic	90	1-11	45	1.7	200 - 5 000
ACQUITY APC XT 125	organic	90	1-11	125	2.5	1 000 - 30 000
ACQUITY APC XT 200	organic	90	1-11	200	2.5	3 000 – 70 000
ACQUITY APC XT 450	organic	90	1-11	450	2.5	20 000 – 400 000
ACQUITY APC XT 900	organic	90	1-11	900	2.5	300 000 - 2 000 000
ACQUITY APC AQ 45	aqueous	45	1-9	45	1.7	200 - 5 000
ACQUITY APC AQ 125	aqueous	45	1-9	125	2.5	1 000 - 30 000
ACQUITY APC AQ 200	aqueous	45	1-9	200	2.5	3 000 – 70 000
ACQUITY APC AQ 450	aqueous	45	1-9	450	2.5	20 000 – 400 000
ACQUITY APC AQ 900	aqueous	45	1-9	900	2.5	300 000 - 2 000 000


Advanced Polymer Chromatograph (APC) / Ultra High Performance Liquid Chromatograph (UHPLC)

1000 bar

**Sample
Manager**

**Diode-array
UV-Vis
detector
(0.5 μ l)**

**Refractive
Index
detector
(1.3 μ l)**




**Thermostat I
for small columns
(4.6 x 150 mm)**

**Thermostat II
for large columns
(7.8 x 300 mm)**


**Quaternary
Solvent
Manager**

Waters Empower 3 software

Behind doors...




Quaternary Solvent Manager




1. Vacuum degasser, four chambers
2. Proportioning valve
3. Pump I
4. Pump II
5. Mixer

Behind doors...




Sample Manager with flow-through needle



Two plates for 96 samples in 2-mL vial holders
Temperature range: 4-40 °C

Behind doors...



Thermostat I for small columns (4.6 x 150 mm)




- 1. column: 450 Å
- 2. column: 200 Å
- 3. column: 125 Å


0.2 μm Filter

Temperature range: 4 – 90 °C


Lignocellulose structure




Lignin structure, most common bonds (β -O-4, β - β) and three types of monomer groups (S, G, H)



4-hydroxyphenyl (H) group




Guaiacyl (G) group




Syringyl (S) group

Lignin valorization


Non-catalytic thermochemical depolymerization




Catalytic depolymerization



Lignin types and conditions of extraction / production




APC chromatogram and UV spectrum of Sodium lignosulfonate (Aldrich)



3 APC QT columns:
450Å + 200Å + 125Å
T = 30°C;
0.5 ml/min
0.3M NaNO₃/MeOH (80:20)


Sample: 1000 µg/mL lignosulphonate; at 280 nm, Sampling rate: 20 pt/s (Hz); Injection volume: 50 µL;

Molecular weight distribution of LB lignin and sequential solvent fractions determined by APC



	LB	EtOAc	EtOH	MeOH	Acetone	Residue
Mn:	700	300	500	800	1300	2600
Mw:	5400	400	900	2500	4300	20500
Mz:	34500	700	2200	6100	8800	62100
D:	7.85	1.29	1.97	2.94	3.38	7.97

APC elution patterns of guaiacol, β -O-4 dimer and guaiacol-based oligomer model



APC XT columns:
125 Å and 45 Å
85 °C;
0.5 ml/min
DMSO + 0.5 % LiBr;
UV at 280 nm

Thank for

Barthos Róbert
Lónyi Ferenc
Novodárszki Gyula
Solt Hanna
Szabó Blanka

Szegedi Ágnes
Vikár Anna
Valyon József
Tuba Róbert

Wellischné Farkas Ágnes
Fekete Miklós

Thank you for your kind attention!

Acknowledgement

European Regional Development Fund (Interreg, SKHU/1902/4.1/001/Bioeconomy)
www.skhu.eu



Building Partnership

